

PREVENTA SCHAMPI PILLOW

Give your bed an additional sitting function



Increased well-being and self-esteem

Eases in care situations for staff

Can prevent decubitus prophylaxis

Can give relief from respiratory distress

” *An easy-to-use aid that provides increased well-being and self-esteem for patients who have difficulty sitting up by them self.*

In places such as homes for the elderly it is quite a common occurrence that the person wants to sit on the edge of the bed, but is unable to. The Schampi pillow has been designed to solve this issue. The pillow is placed behind the back of the person, once he or she has been helped into a sitting position. The pillow is prevented from sliding backward by securing the side rail of the bed behind the person. It provides a secure and comfortable positioning thanks to the contoured back rest.

The user's self-esteem

Being able to sit on the edge of the bed is important for an individual's self-esteem. Becoming a "prisoner" in the supine position can be debilitating and depressing. It can provide increased well-being and feel stimulating to sit on the edge of the bed and eat or read the paper without assistance.

Moreover, many patients who suffer from anxiety/respiratory distress experience relief when they can sit on the edge of the bed for a moment. This could be applicable to care recipients with lung disease, such as asthma and COPD, or heart failure.

The Schampi pillow also has positive effects on pressure ulcer prevention, as it increases positioning options. Variation is made possible/encouraged. One-sided body positioning (supine) is avoided.

Provides relief for care staff

The Schampi pillow has been developed with the patient in mind, but also for the convenience of the care staff by releasing resources and improving ergonomics. The Schampi pillow provides support that can be helpful in feeding, dressing/undressing and other care duties.

Areas of use

All types of hospitals, homes for the elderly, service apartments, special needs housing, palliative care units and rehabilitation units. Other areas such as care homes for disabled and handicapped individuals, community-based housing, home care, and maternity wards.



Schampi pillow with armrest



Schampi pillow Small



Schampi pillow without armrest

Article no.	Description	Size (LxWxH)	Color
13-000100-001	Low back with armrest, PU cover	85x70x40 cm	blue/beige
13-000100-002	Smaller width with high back, PU cover	85x58x53 cm	blue/beige
13-000000-001	Low back without armrest, PU cover	40x50x38 cm	blue/beige



Cover

The pillow is available with a cover made of PU material. The cover is soft, comfortable and easy to clean. The cover is easy to remove for cleaning purposes (Velcro tape, see picture below).

Recycling

The Schampi pillows are fully recyclable and should be sorted as combustible waste. No dismantling is required.



Please read the instruction manual before using the product.

