# OPTICELL AIR1S

Hybrid mattress with pump for pressure ulcer category 4 TruePressure high comfort and extra low pressure

TrueAlternating ensures a smoother air exchange

OptiTex<sup>™</sup> cover with excellent features

*An alternating hybrid mattress with pump, integrated air cells, and hygiene cover.* 

OPTICELL is Järven Health Cares product series of pressure-relieving air mattresses.

## OPTICELL AIRIS

## AirIS with TruePressure is an alternating hybrid mattress with pump OC20 V2 designed to meet the stringent demands of healthcare and the needs of the patient.

Its interconnected air cells distribute the air pressure over the mattress providing both high comfort and extra low back pressure\*. We call this **TruePressure**.

The mattress has been designed using carefully selected components to provide the conditions for immersion and envelopment. Two factors that research has shown to be of great importance for effective and even pressure relief.\*\*

#### Features – Pump and cells

COMFORT FUNCTION with six different comfort settings. TRANSPORT MODE – The mattress can be set in transport mode where it maintains maximum air pressure for at least 12 hours ALARM FUNCTION – The air pump has an alarm function with both light and sound alarms that give warnings for low pressure in the mattress. HEEL RELIEF – The mattress is fitted with a relief system with an integrated airlock that enables the easy disconnection of cells. EXTRA LONG POWER CABLE – We realise that the distance to the nearest power socket may vary which is why the mattress comes with a 6 metre long cable.

#### OptiTex<sup>™</sup> cover

**GOD HYGIENE** – The OptiTex<sup>™</sup> cover can be cleaned with disinfectant and can also be machine washed at 95°C. The cells and hoses can be rinsed clean with warm water (max 70°C) or with disinfectant.

**DRIP PROTECTION WATERPROOF ZIP** on four sides. All seams are welded. **THE UNDERSIDE HAS ANTI-GLIDE PROTECTION** with four fixing strips. **FOUR HANDLES** to facilitate handling of the mattress.

#### Foam configuration

**WATERLILY FOAM** – A high-quality foam with a low environmental impact. **HYPERSOFT FOAM** – A really soft foam that moulds to the body. Together, these foams facilitate immersion and envelopment while achieving a high level of comfort.

	Pressure values, Air1 (without pump) hybrid mattress				
Raised mattress	Large impression body 300 N	Large impression body 600 N	Small impression body 40 N		
20 cm	26	45	42		
18 cm	26	47	44		

Reference values						
Impression body	Extra low	Low	Medium	High		
Large impression body 300 N	< 30.0	30.0-34.9	35.0-39.9	40.0-44.9		
Large impression body 600 N	< 45.0	45.0-52.9	53.0-60.9	61.0-68.9		
Small impression body 40 N	< 45.0	45.0-52.9	53.0-60.9	61.0-68.9		

\* RISE Research Institutes of Sweden O100560-1129023

\*\* What is new in our understanding of pressure injuries, Amit Gefen and Esa Soppi 2020.

\*\*\* Pressure ulcer categories defined according to the International Pressure Ulcer Classification system. Guidelines produced by the European Pressure Ulcer Advisory Panel (EPUAP), National Pressure Ulcer Advisory Panel (NPUAP) och Pan Pacific Pressure Injury Alliance (PPPIA).



### Product**specifications**

**STANDARDSIZES** 200/210/220 x 80/85/90/105/120 cm. OptiCell Air1 comes in different mattress heights and maximum loads. See table with data for each model. For other sizes please contact Järven Health Care.

#### PRESSURE ULCER CATEGORY 1-4\*\*\*

 $\mbox{PUMP OC20 2.0}$  is an advanced pump with sensors that adapts to the user's weight Sound level: 21 db. 230 V / 50 Hz.

**CELLS** Our unique replaceable SweCells™ are made of TPU (Thermoplastic Polyurethane) and are softer than traditional PU cells. This means that they mould more easily to the body with less rustling noise.

COVER of bi-elastic double-coated textile with soft elastic PU.

WARRANTY 2 years for manufacturing defects.

FIRE STANDARD SS 876 00 01, EN 597-1 & 2.

**PATIENT WEIGHT** OptiCell Air1 comes in different mattress heights and maximum loads. See table with data for each model.

#### **RECYCLING** See User Manual

 WEIGHT with pump
 11.8 kg (200x80x16 cm)

 13.2 kg (200x80x16 cm)

LABELLING year of manufacture, month, week, EAN, care instructions, standard compliance and placement on bed.





Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN Phone: + 46 (0)660 - 22 28 50 | Fax: + 46 (0)660 - 824 10 www.jarven.se | healthcare@jarven.se