



FLEX PAD

Our pressure ulcer prevention heel pad

The pad is used with patients who have lower limb pressure-sensitive areas. It can be used at the first signs of a pressure sore to prevent further injury and reduce the risk of exacerbating the pressure sore. The Flex Pad is also suitable for patients who have been given spinal anaesthesia or who suffer from pain problems.

- **TWO VERSIONS.** The pad is available in a choice of two models, one for heavier patients (85-130 kg) with a bevelled edge and one for lighter patients with a plane edge.
- **UNIQUE SURFACE LAYER LENTEX®.** The Lentex® surface layer is laminated directly onto the foam core, which means that it is seamless and zipperless. The surface layer is white to make it easier to detect contaminants.

PRODUCT SPECIFICATIONS

Size	Plane edge: 80 × 50 × 4 cm Bevelled edge: 80 x 55 x 7 cm or 80 x 55 x 12 cm
Materials	Carefully selected viscoelastic foam. Cover of the Lentex®.
Guarantee	12 months for manufacturing defects. IMPORTANT The guarantee does not cover damage caused by physical impact with a pen or other sharp object.
Fire Safety	SS 876 00 01, test method NT-Fire 037 A+B.
Recycling	The Preventa Flex Pad is fully recyclable and should be sorted as combustible waste. No dismantling necessary.

