

# MELISSA BASIC

## Reversible and rotatable care mattress with optional hygiene cover

With Melissa Basic you get an affordable and good basic mattress with good durability and maintained quality. Melissa Basic is designed to help prevent pressure ulcers up to and including category 1.

- **HIGH RESILIENT COLD FOAM.** Melissa Basic has a high resilient cold foam that provides a durable and stable healthcare mattress.
- **WATERPROOF COVER IN TWO VARIANTS.** Melissa Basic comes with a noiseless RX cover. The cover can be supplied with a handle or evacuation function.
- **ZIPPER ON 3 OR 4 SIDES.** Depending on the choice of cover, the cover has zippers on three or four sides. All covers have a drip shield over the zippers to facilitate hygiene work.



# MELISSA BASIC

## AN EXCELLENT STANDARD MATTRESS WITH A REVERSIBLE AND ROTATABLE FOAM CORE.


The mattress has a homogeneous foam core which allows it to be flipped and turned. Something that allows a more even load over its surface. The Melissa Basic is designed to be a quiet worker over many years.

### Waterproof cover with good properties

Melissa Basic comes with RX Cover. A hygienic cover made of 2-stretch polyurethane fabric that is breathable and noiseless. The cover is available in three models:

- RX (supplied as standard).
- RXH with handles that make it easier to manage the mattress.
- RXI with built-in evacuation function.

Covers have a zipper with a drip shield on 3 sides except for cover RXI which zips with a drip shield on 4 sides.



### STANDARD FOR THE OBJECTIVE EVALUATION OF THE PRESSURE-RELIEVING PROPERTIES OF MATTRESSES.

In order to obtain an objective assessment of the pressure relieving properties of mattresses, RISE has developed a standard where the actual counterpressure of a mattress is measured by means of an indentational body. The measurement method is described in the Swedish standard SS 876 00 13.

Melissa Basic has been tested by RISE\* according to method SS 876 00 13 and is pressure relieving up to wound and including category 1.

**The test results are presented below.**

Load capacity	Extra low	Low	Medium	High
Large indenter 300 N	< 30,0	30,0 - 34,9	35,0 - 39,9	40,0 - 44,9
Large indenter 600 N	< 45,0	45,0 - 52,9	53,0 - 60,9	61,0 - 68,9
Small indenter 40 N	< 45,0	45,0 - 52,9	53,0 - 60,9	61,0 - 68,9

Mattress	Large indenter 300 N	Large indenter 600 N	Small indenter 40 N
Melissa Basic	38	52	68

\*RISE (Research Institutes of Sweden), a Swedish independent research institute - [www.ri.se](http://www.ri.se)

## PRODUCT SPECIFICATIONS



<b>Materials</b>	12 cm cold foam 35 kg/m <sup>3</sup>
<b>Sizes</b>	200x80/85/90/105/120x12 cm (For other sizes contact Järven Health Care.)
<b>Wound category</b>	Up to and including category 1
<b>Mattress weight</b>	approx. 7 kg depending on length and width

<b>User weight</b>	Up to and including 180 kg.
<b>Fire safety</b>	SS 87600 01, normal risk (test method EN 597-1 & 2)
<b>Recycling</b>	Mattress with cover is fully recyclable. Sort as combustible waste. Disassembly is not required.
<b>Marking</b>	Year of manufacture, month, day, EAN, care instructions, contact information, standards met and bed location.



Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN  
Tel: +46 (0)660 - 22 28 50 | Fax: +46 (0)660 - 824 10  
[www.jarven.se](http://www.jarven.se) | [healthcare@jarven.se](mailto:healthcare@jarven.se)