



MELISSA PARAD- *Mattress overlay*

A high-quality mattress overlay developed to create safe healthcare environments

Experience high comfort and prevent pressure ulcers with the Melissa Parad mattress overlay. A 6 cm overlay mattress with pressure relieving properties, this mattress is a good complement to your existing homecare mattress. In addition, with the ability to flip and rotate the core, you get long-term use while maintaining quality.

- **HIGHLY RESILIENT HYPERSOFT FOAM.**
The Melissa Parad mattress overlay is made from high resilient hyper-soft foam which provides excellent pressure relief.
- **UP TO AND INCLUDING WOUND CATEGORY 3.**
Together with a healthcare mattress, the Melissa Parad mattress overlay is designed to help treat patients with pressure ulcers up to and including wound category 3.
- **WATERPROOF COVER WITH ELASTIC.**
The mattress overlay comes with a noiseless RX cover. The cover is supplied with elastic bands at the corners to keep the mattress overlay in place.
- **REVERSIBLE FOAM CORE.**
The foam core of this mattress overlay can be both flipped and rotated without altering its properties. By doing so, the lifespan of your Melissa Parad mattress overlay will be extended.




MELISSA PARAD - *Mattress overlay*

A GREAT COMPLEMENT TO YOUR HEALTHCARE MATTRESS.

Melissa Parad mattress overlay is designed to provide pressure relief and comfort. The pressure relieving mattress provides relief for sensitive areas such as shoulders, heels and elbows, eliminating the need for a separate heel zone. The two-stretch hygienic cover is also designed to reduce the risk of shearing.

Waterproof cover with good properties

Melissa Parad mattress overlay comes with a cover from the RX series. A series of covers made of 2-stretch polyurethane fabric that is both breathable and noiseless. Depending on the choice of cover, the cover has zippers on three or four sides. All covers have a drip shield covering the zipper to facilitate hygiene.



STANDARD FOR THE OBJECTIVE EVALUATION OF THE PRESSURE-RELIEVING PROPERTIES OF MATTRESSES.

In order to obtain an objective assessment of the pressure relieving properties of mattresses, RISE has developed a standard where the actual counterpressure of a mattress is measured by means of an indentational body. The measurement method is described in the Swedish standard SS 876 00 13.

Melissa Parad mattress overlay has been tested by RISE* according to method SS 876 00 13 and is pressure relieving up to wound category 3.

The test results are presented below.

Load capacity	Extra low	Low	Medium	High
Large indenter 300 N	< 30,0	30,0 - 34,9	35,0 - 39,9	40,0 - 44,9
Large indenter 600 N	< 45,0	45,0 - 52,9	53,0 - 60,9	61,0 - 68,9
Small indenter 40 N	< 45,0	45,0 - 52,9	53,0 - 60,9	61,0 - 68,9

Mattress	Large indenter 300 N	Large indenter 600 N	Small indenter 40 N
Melissa Parad, mattress overlay	24	43	42

*RISE (Research Institutes of Sweden), a Swedish independent research institute - www.ri.se

PRODUCT SPECIFICATIONS



Materials	6 cm cold foam 50 kg/m ³	User weight	Up to and including 180 kg.
Sizes	200/205/210 x 80/85/90/105/120 x 6 cm (For other sizes contact Järven Health Care.)	Fire safety	SS 87600 01, normal risk (test method EN 597-1 & 2)
Wound category	Up to and including category 3.	Recycling	Mattress with cover is fully recyclable. Sort as combustible waste. Disassembly is not required.
Mattress weight	approx. 5 kg depending on length and width.	Marking	Year of manufacture, month, day, EAN, care instructions, contact information, standards met and bed location.



Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN
Tel: +46 (0)660 - 22 28 50 | Fax: +46 (0)660 - 824 10
www.jarven.se | healthcare@jarven.se