


OPTICELL GENTLE



Up to
and including
category 3

Auto-adapts
to the patients
weight

Deflates rapidly
for emergency
CPR use

OptiTex™ cover
with innovative
features

*A pain relieving static air mattress with hygiene cover,
automatic pump and interchangeable cells.*

OPTICELL GENTLE

The OptiCell Gentle is an 18cm, constant low pressure (CLP), full replacement air mattress system with an integrated foam core base which can be used as an aid in the prevention and treatment of pressure ulcers up to and including category 3.

This mattress system includes a fully stretchable and removable hygiene cover with welded seams and drip protection. It also includes 20 individual air cells which can be easily removed.

The OptiCell Gentle includes ventilated air cells which provide a user-friendly microclimate which reduce the risk of skin maceration while keeping the top cover dry.

An easy to use, energy efficient and silent OC31 pump, is also included. Inflation time is approx. 15 minutes.

FEATURES – Pump and mattress cells

CPR VALVE – Ensures rapid deflation in an emergency situation requiring CPR.

ALARM FUNCTION – Uses both light and sound alarms to warn of low pressure in the mattress or of power outages.

HEEL RELIEF – Equipped with SweCells™ and an air release system which enables easy deflation of individual cells by means of an integrated air-lock, to completely eliminate the pressure on heels.

EXTRA LONG POWER CABLE – Comes with a 6-metre cable to ensure suitable access to the nearest power outlet.

FEATURES – OptiTex™ cover

INTEGRATED WATERPROOF CABLE HOLDER to prevent bacteria spreading, and keep the electrical cord out of the way and protects the air hoses.

HYGIENIC – The entire cover (both upper and lower parts) can be cleaned with disinfectants or machine washed at maximum 95° C.

Cells and hoses are easily removed and rinsed clean with disinfectants or with warm water at maximum 70° C.

LOW AIR LOSS FUNKTION tillåter luft att komma ut från cellerna samtidigt som det är vattentätt. Detta för att minska risken för hudmaceration och andra komplikationer för brukaren.

FLUID PROTECTION FLAP AND WATERPROOF ZIP on all four sides, with welded seams.

ANTI-SLIP UNDERSIDE grips to bed base in addition to four fastening straps.

FOUR TRANSPORTATION HANDLES TO ENABLE EASY MOVING

Produktspecifikationer

STANDARD SIZES 200/210/220 x 80/85/90/105/120 X 18 cm.

For other sizes please contact Järven Health Care.

PRESSURE ULCER CATEGORY 1–3*

CELLS The unique interchangeable SweCells™ are made of thermoplastic polyurethane (TPU) and are softer than traditional polyurethane (PU) cells. This means they conform to the patient's body more effectively, giving increased side support, with less noise.

COVER Multi-stretch, waterproof, breathable and washable OptiTex™ bi-elastic double coated textile with soft elastic PU.

PUMP OC31 is an advanced pump with sensors that adapt to the patient's weight. A dual-function indicator lamp shows when the pump and mattress are ready to use, or if there is a fault requiring service.

Sound level: 21 db. 230 V / 50 Hz.

WARRANTY 2 years on manufacturing defects.

FIRE SAFETY STANDARDS Complies with SS 876 00 01, EN 597-1 & 2.

USER WEIGHT Up to 250 kg.

RECYCLING See instructions.

WEIGHT Mattress 11 kg (200 x 105 cm). Pump 34 kg

LABELLING Date of manufacture, EAN, care instructions, fulfilled standards and placement in bed.
The mattress and cells are PVC-free



* Categories defined by the International Pressure Ulcer Classification system. Guidelines produced by the European US National Pressure Ulcer Advisory Panels (EPUAP and NPUAP), together with the Pan Pacific Pressure Injury Alliance (PPPIA).



OptiTex™ – hygienöverdrag med goda egenskaper



ÖVERDRAGETS
VATTENTÄLIGHET
ISO 811-1981



LUFTGENOM-
SLÄPPLIGHET
ISO 11092



ÖVERDRAGETS
BRANDSTANDARD
BS EN 597-1 & 2



ÖVERDRAGETS
STRÄCKTÄLIGHET
EN ISO 13934-1 1999



HYDROSTATISK
FÖRMÅGA
ISO 811-1981



Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN
Telefon: 0660 - 22 28 50 | Telefax: 0660 - 824 10
www.jarven.se | healthcare@jarven.se